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 Packaging unit

Purpose:

The **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tubes with Cuff are intended for patients with swallowing disorders. They feature a cuff for blocking and a universal adapter that can be used for connection to a speech valve or to accommodate a paediatric HME (heat-moisture exchanger). The universal adapter is extended towards the front depending on the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube version. The **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tubes are also available with an extension of the universal adapter to simplify connections to it in difficult anatomical situations.

Description:

Before the tube is inserted, the cuff must be inflated to make sure it does not leak using a cuff pressure measuring device. If the cuff is leaktight, allow the air to escape; the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube can then be inserted by a physician or nursing staff. To ensure easy insertion of the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube into the tracheostoma it should be treated with a suitable lubricant. When blocking the cuff, make sure not to use pressure greater than 25 mmHg (equivalent to 34 cmH₂O) over a longer period. That is why the cuff pressure should be checked at regular intervals using a cuff pressure measuring device. If the tube becomes clogged due to pronounced secretion, the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube must be removed and cleaned. Then use a tracheal suction catheter to suction-clean it through the outer tube. The **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube is then reinserted.

PRIMA-STOM II PED C/LC Paediatric Tracheal Tubes are one-patient products. They are made of medically certified plastic materials. The tube adapts to the anatomical features of the trachea due to the body heat. The maximum period of use is 29 days. After a maximum of 6 days in place, the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube should be removed so that the trachea can be examined. This makes possible early recognition or avoidance of stenoses, for example. Before reinserting the **PRIMA-STOM II PED C/LC** Paediatric Tracheal Tube it must be cleaned as explained in the section "Cleaning".

Indication:

PRIMA-STOM II PED C/LC Paediatric Tracheal Tubes are particularly suited for use in children suffering from swallowing disorders (dysphagias) of all aetiologies and levels of severity to protect the lower respiratory tract. This patient group also includes ventilated children during weaning (that is while they are adapting during the post-ventilation period) and following tumour removal or radiation treatment in the head and neck region.

Cleaning:

The tube is cleaned without chemical additives in a saline solution at a maximum of 65°C. Following this we recommend disinfection of the internal and external tubes using a suitable disinfectant following the instructions for use. Following cleaning and disinfection, the cuff must be tested to make sure it does not leak.

Safety information:

The product must not be used if the packaging is damaged. Contamination with candida is possible if the cleaning instructions are not followed, the internal tube could then no longer have a snug fit in the external tube. Excessively high blocking pressure in the cuff may cause pressure necroses in the trachea.

Completely unblock the maintenance cannula before introducing it into the trachea. Prior to blocking, the maintenance cannula should not be retracted to prevent obstruction of the inner diameter by the block.

Caution: Leaving the protective cap on the pilot balloon during use may result in a loss of cuff pressure. The protective cap can be discarded once the tube has been unpacked. It is only used for initial sterilisation by the manufacturer.

The product is not suitable for use with laser devices.

Serious incidents involving the maintenance cannulas must be reported to the manufacturer and the competent authority.

Sterility Information:

The products have been sterilised using ethylene oxide. Sterility is guaranteed only if the packaging is unopened and free of damage or tears. The products may not be re-sterilised.

Product disposal:

Disposal of products must be carried out in compliance with current national regulations applying to waste products.

Caution: Warranty, guarantee and liability claims against the manufacturer will be nullified if any changes are made to the product by users or if the product is used outside of the scope of applications described in these instructions.